

# VersaGlue<sup>®</sup>PX

## Laboratory Bonding System

dental**lab**products<sup>®</sup>

Best of the Best  
PRODUCTS

2008

### Performance Extended (PX) Dental Industry Adhesives

#### **VERSATILE**

Bonds all dental materials, including plaster, metal, acrylic and porcelain, wet or dry.

#### **RELIABLE PRECISION**

Precision applicator allows for precise material placement with minimal waste.

#### **RAPID STRUCTURAL BOND**

Highly refined and pure Ethyl formula offers optimal tensile strength resulting in instant, powerful, and impact-free bonding.

#### **EXTENDED SHELF-LIFE**

Up to 3 years unopened and refrigerated.  
Up to 6 months opened.

Assembled in the USA of  
foreign and US components



**HarvestDental**  
LABORATORY PRODUCTS

Toll-Free: 800-706-7599    [www.harvestdental.com](http://www.harvestdental.com)



## VersaGlue®PX-4 LOW VISCOSITY

### Properties:

Low viscosity, thin liquid flows easily into tight and hard to reach areas.

Set time is 1-5 seconds.

### Applications:

Gluing articulating systems, sealing dies and hardening models, repairing stone dies, broken denture teeth and popped veneers.

Item 2025, 1 oz

## VersaGlue®PX-40 MEDIUM VISCOSITY

### Properties:

Medium viscosity glue exhibits exceptional gap filling properties.

Set time is 10-15 seconds.

### Applications:

Dowel pin placement, fixing broken models & dies, assembling crowns & bridges for soldering.

Item 2026, 1 oz

## VersaGlue®PX-400 HIGH VISCOSITY

### Properties:

High viscosity glue for closing large gaps. Thick and strong for high stress areas, allows more time for positioning.

Set time is 30-60 seconds.

### Applications

Dowel pin placement, Filling undercuts & voids, stabilizing abutment segments.

Item 2027, 1 oz

**FOR PROFESSIONAL DENTAL USE ONLY.  
NOT FOR INTRA-ORAL USE.**

**KEEP OUT OF REACH OF CHILDREN**



**HarvestDental**  
LABORATORY PRODUCTS

Assembled in the USA of foreign and US components.

Toll-Free: 800-706-7599  
www.harvestdental.com